- ⌚ Thread a PVA bag (Ref. W56) and inject the resin (Ref: W0905)
- ⌚ Before the end of stratification tighten with a link in the throat of the anchor
- ⌚ After taking full resin, remove the dough in the anchor and any excess resin

With tubular carbon stockinette W57GA

- ⌚ Isolate the positive with a PVA bag (Ref: W56 or W56R)
- ⌚ Cover with a first polyamide fabric (Ref: W52 or W53)
- ⌚ Cut several slices of tissue carbon (ref: W57TA or W57SE) and position on the future site of the anchor to cement
- ⌚ We recommend using repositionable adhesive (ref: Z35) for this operation
- ⌚ Position the anchor and check the alignment
- ⌚ Cut a length of carbon stockinette (ref: W57GA) and put on the positive
- ⌚ Use a link (with Kevlar type) and tighten the braid tube in the throat of the anchor to cement
- ⌚ Trim excess braid or flatten it on the positive
- ⌚ Cover with a Polyamide tubular stockinette (Ref: W52 or W53) or tubular sheath for finishing (Ref: W54A-25)
- ⌚ Thread a PVA bag (ref: W56) and inject the resin (ref: W0905)
- ⌚ Before the end of stratification tighten with a link in the throat of the anchor
- ⌚ After taking full resin, remove the plasticine and any excess resin

Adjustments

Use screws of the adapter (see technical instructions for the adapter).

Restriction on use / Precautions for use

- ⌚ The adaptor must be used within the recommended weight limits and in accordance with the assessment made by an approved ortho-prosthetist depending on the patient concerned
- ⌚ Product maintenance is subject to an assessment made by an approved ortho-prosthetist depending on the conditions of use
- ⌚ All assembly and adjustment work must be carried out by qualified personnel

Cleaning / Storage

- ⌚ Store the product in a dry place, in its original packing (on which the batch number is marked)
- ⌚ The frequency of inspections to check the equipment is subject to the assessment made by the ortho-prosthetist depending on the patient concerned
- ⌚ Use a damp cloth to clean thea product

Guarantee

The guarantee is valid for one year as from the date of purchase of the product, if the product has been used in accordance with the manufacturer's instructions.

N.B. In the event of a complaint, please enclose a return form correctly filled in; incomplete forms cannot be dealt with.

EC markings

Groupe G2M hereby declares, in its quality as the manufacturer, that all its parts comply with directive 93/42/EC.

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Designation

Titanium tibial anchor for cementing, with 4 prongs with rotating receiver adapter.

- Weight : 101 g

Purpose

Structural element used to make endoskeletal prostheses of the lower limbs for patients having average levels of activity and whose weight does not exceed 136 kg.

Compatibilities

The tibial anchor for cementing is compatible with all the adaptors in the titanium range.

Fitting

Instructions for stratification

- Preparation of the anchor CA115Ti
- Remove the 3 screws of the anchor, remove the adapter
- Fill the interstices with plasticine (ref. ZOP516)

Using W57SE or W57TA carbon tissue

- Isolate the positive with a PVA bag (Ref: W56 or W56R)
- Cover with a first polyamide fabric (Ref: W52 or W53)
- Cut a piece of carbon (ref: W57SE or W57TA) and apply over the entire surface of the positive
- Cut slits in the surplus from the top and back on the positive
- Cut several slices of tissue carbon and position on the future site of the anchor to cement
- Cut and position a piece of carbon cloth to cover the whole positive. We recommend using repositionable adhesive (ref: Z35) for these operations
- Position the anchor and check the alignment
- Cut a strip of cloth to cover the carbon branches of the anchor
- Cut pieces of carbon to strengthen the parts you want (ischium, supports sub-patellar, medial, etc..) Positioning with repositionable adhesive
- Cut a second piece of carbon cloth to cover the entire.
- Use a link (with Kevlar type) and tighten down the throat of the anchor to sink
- Cover with a cloth Polyamide (Ref: W52 or W53) or cloth trim (Ref: W54A-25)