Fitting


Adjusting connector orientation

- Screw the connectors fully onto the CA131STFX/STMX 4 hole plate
- Adjust the rotation by unscrewing the connector ¼ turn at most
- To maintain the desired orientation, tighten the locking screw (3) on the plate

Adjustment

Tilt amplitude : +/- 9° compared to the reference axis.
Set the desired tilt using the 4 screws, in accordance with the limits set.
Lock the assembly using (Z31-243) locking compound on the screws.

Spare parts

		CA133STFX	CA133STMX
1	CHC screw M8 L16	V03I-0816	V03Z-0816
2	FHC screw M6 L16	V05I-0616	V05Z-0616
3	CHC screw M5 L30	V02I-0530	V02FZ-0530

Restriction on use / Precautions for use

Tightening torques (use key (ZOC102))

Screws	Torques (Nm)	Screws	Torques (Nm)	Screws	Torques (Nm)
V03I-0816	13	V05I-0616	9	V02FZ-0530	6
V02I-0530	6	V03Z-0816	15	V05Z-0616	9

Cleaning / Storage

- Store the product in a dry place, in its original packing (on which the batch number is marked)
- The frequency of inspections to check the equipment is subject to the assessment made by the ortho-prosthetist depending on the patient concerned

Guarantee

The guarantee is valid for one year as from the date of purchase of the product, if the product has been used in accordance with the manufacturer's instructions.


N.B. In the event of a complaint, please enclose a return form correctly filled in, incomplete forms cannot be dealt with.

EC markings

Groupe G2M hereby declares, in its quality as the manufacturer, that all its parts comply with directive 93/42/EC.

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